GB

The micro-current beauty instrument brings personal training to a new level of beauty. It provides a series of customized micro-current facial exercise modes. The micro-current contact ball can directly conduct micro-current to the skin to help generate collagen, repair skin elasticity and tighten facial contours.

Material: ABS plastic Dimensions: 57x52x24 mm Weight: 35 grams

Battery: 3.7V at 180mA Frequency: 185 Hz Noise level: less than 50 decibels

Personal trainer for facial exercises

Beauty double head instrument function

1. Short press the start button, the beauty function will be turned on, micro current level 1, and then short press micro current level 2, a total of 3 cycles.
2. Double click vibration to turn it on, and then double click vibration to turn it off.
3. Press and hold for 2 seconds to turn off the instrument.
4. When charging, the red light is on and the full red light is off.

Combining precise micro-current and vertical skin-penetrating sonic pulsation technology, facial muscle cells can be activated in just a few minutes. The bottom of the muscle can generate low-voltage currents from bionic micro-currents, mimicking the body's own bio-electricity, and reaching the bottom of the muscles gently.

Deeply activate muscle movement and stimulate firming energy.

Make your skin firmer and have a more radiant complexion-rejuvenate your skin.

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